

**UNCLASSIFIED**

---

**AD 295 672**

*Reproduced  
by the*

**ARMED SERVICES TECHNICAL INFORMATION AGENCY  
ARLINGTON HALL STATION  
ARLINGTON 12, VIRGINIA**



---

**UNCLASSIFIED**

NOTICE: When government or other drawings, specifications or other data are used for any purpose other than in connection with a definitely related government procurement operation, the U. S. Government thereby incurs no responsibility, nor any obligation whatsoever; and the fact that the Government may have formulated, furnished, or in any way supplied the said drawings, specifications, or other data is not to be regarded by implication or otherwise as in any manner licensing the holder or any other person or corporation, or conveying any rights or permission to manufacture, use or sell any patented invention that may in any way be related thereto.

SSD-TDR-62-203

**295 672**


**AN ACETYLATION METHOD FOR  
DETERMINING MIXTURES OF EITHER  
HYDRAZINE-1,1-DIMETHYLHYDRAZINE  
OR  
MONOMETHYLHYDRAZINE-1,1-DIMETHYLHYDRAZINE**

TECHNICAL DOCUMENTARY REPORT NO. SSD-TDR-62-203

DECEMBER 1962

ROCKET RESEARCH LABORATORIES  
CHEMICAL AND MATERIALS BRANCH  
EDWARDS, CALIFORNIA

HEADQUARTERS  
SPACE SYSTEMS DIVISION  
AIR FORCE SYSTEMS COMMAND  
UNITED STATES AIR FORCE



## NOTICES

When Government drawings, specifications, or other data are used for any purpose other than in connection with a definitely related Government procurement operation, the United States Government thereby incurs no responsibility nor any obligation whatsoever; and the fact that the Government may have formulated, furnished, or in any way supplied the said drawings, specifications, or other data, is not to be regarded by implication or otherwise as in any manner licensing the holder or any other person or corporation, or conveying any rights or permission to manufacture, use, or sell any patented invention that may in any way be related thereto.

Qualified requesters may obtain copies of this report from the Armed Services Technical Information Agency, (ASTIA), Arlington Hall Station, Arlington 12, Virginia.

SSD-TDR-62-203

AN ACETYLATION METHOD FOR DETERMINING  
MIXTURES OF EITHER HYDRAZINE-1, 1-DIMETHYLHYDRAZINE  
OR MONOMETHYLHYDRAZINE-1, 1-DIMETHYLHYDRAZINE

TECHNICAL DOCUMENTARY REPORT  
NO. SSD-TDR-62-203


DECEMBER 1962

by

HUGH E. MALONE

ROBERT A. BIGGERS

This report has been reviewed and approved.



HAROLD W. NORTON  
Colonel, USAF  
Commander

Rocket Research Laboratories  
Chemical and Materials Branch  
Edwards, California

### ABSTRACT

A method is presented for determining mixtures of either hydrazine-1, 1-dimethylhydrazine ( $\text{N}_2\text{H}_4$ /UDMH) or monomethylhydrazine-1, 1-dimethylhydrazine (MMH/UDMH). This method is based on the difference in reaction rates for acetylating hydrazines. Both  $\text{N}_2\text{H}_4$  and MMH react immediately while the acetylation of UDMH proceeds slowly. The reactions are conducted in acetonitrile-acetic acid medium. Two titrations are required: one measures the total basicity of the hydrazine mixtures; the other titration measures UDMH after the  $\text{N}_2\text{H}_4$  or MMH has reacted with  $\text{Ac}_2\text{O}$ .

## TABLE OF CONTENTS

	<u>Page</u>
INTRODUCTION . . . . .	1
BACKGROUND . . . . .	1
EXPERIMENTAL . . . . .	2
Effect of Acetic Anhydride . . . . .	3
Effect of Water . . . . .	4
Potentiometric and Colorimetric Endpoints . . . . .	5
PROCEDURE . . . . .	6
Preparation of the Sample . . . . .	6
Determination of Total Basicity of the Hydrazine- UDMH Mixture . . . . .	6
Determination of 1, 1-Dimethylhydrazine . . . . .	6
CALCULATIONS . . . . .	6
REAGENTS . . . . .	7
RESULTS . . . . .	7
SUMMARY . . . . .	9
BIBLIOGRAPHY . . . . .	14

## LIST OF ILLUSTRATIONS

### Figure

I.	Effect of Time and Temperature on the Reaction Rate of UDMH with Acetic Anhydride . . . . .	11
II.	The Effect of Acetic Anhydride on Hydrazines . . . . .	12
III.	Typical Titration Curve for Both Titrations . . . . .	13

# AN ACETYLATION METHOD FOR DETERMINING MIXTURES OF EITHER HYDRAZINE-1, 1-DIMETHYLHYDRAZINE OR MONOMETHYLHYDRAZINE-1, 1-DIMETHYLHYDRAZINE

## INTRODUCTION.

Hydrazine ( $N_2H_4$ ), monomethylhydrazine (MMH), and 1, 1-dimethylhydrazine (UDMH) are important both individually and in mixtures as rocket fuels. A 50-50 percent mixture of  $N_2H_4$ /UDMH is used as the fuel blend for the Titan II. Though mixtures of MMH/UDMH are discussed herein, they are not presently being used in any weapon system.

The purpose of this report is to present a method for determining mixtures of either  $N_2H_4$ /UDMH and MMH/UDMH based on observed differences in reaction times with acetic anhydride.

## BACKGROUND.

Several methods have been developed for analysis of mixtures of hydrazine-1, 1-dimethylhydrazine. These methods use gas chromatographic techniques (1, 5), standard iodate titration in an aqueous medium (9), or the selective reaction of hydrazine with salicylaldehyde and subsequent titration of basicity with perchloric acid in non-aqueous medium (6). This latter method is presently used in Specification MIL-P-27402. It has the end point disadvantage encountered in most weak base-strong acid titrations.

Acetic anhydride used herein for acetylating hydrazines was used primarily as a solvent for titrations of amides (10), weak bases (3), and nitrogen bases (8). It has been used also to improve the end point in non-aqueous titrations (10). Fritz and Schenk (4) used acetic anhydride for acid catalyzed acetylation of hydroxyl compounds in ethyl acetate.

Although the reactions for acetylating hydrazines are known (2) and reported (7), the authors are unaware of any efforts to determine hydrazine and substituted hydrazines by this method.




## EXPERIMENTAL.

Reactions of hydrazine, monomethylhydrazine and 1, 1-dimethylhydrazine with acetic anhydride. Reaction products are shown in the following equations:


With  $N_2H_4$ -1, 2 diacetylhydrazine:

Equation (1)




With MMH-Monomethyl 1, 2 diacetylhydrazine:

Equation (2)



With UDMH-1, 1-Dimethyl 2 acetylhydrazine:

Equation (3)



The rates of the reaction differ greatly as both the  $N_2H_4$  and MMH are rendered neutral to perchloric acid titration almost immediately upon UDMH

reacts slowly in acetic acid medium. The product of equation (3) occurs only after a considerable time lapse at room temperature. Figure I shows the effect of time and temperature on the UDMH-acetic anhydride reaction. The reaction, as indicated by a decrease in milliliters of  $\text{HClO}_4$  is slow until the temperature is increased.

#### Effect of Acetic Anhydride.

As the quantity of acetic anhydride is increased in samples of  $\text{N}_2\text{H}_4$ , MMH, and UDMH, the milliliters of  $\text{HClO}_4$  titrant required for neutralization decreases much faster with  $\text{N}_2\text{H}_4$  and MMH. The UDMH is slow to react even when a considerable excess of acetic anhydride is added. This effect is shown in Figure II. These curves were plotted from titration, with perchloric acid, of individual samples containing a 1 ml aliquot of each hydrazine compound per 50 ml acetic acid. Duplicate samples were treated with 1 ml aliquots of acetic anhydride samples containing 1.0 ml - 6.0 ml acetic anhydride in 50 ml acetic acid. The two curves for both MMH and  $\text{N}_2\text{H}_4$  substantiate the molar ratio presented in equation (1) and (2).

An analysis, after acetylation, of the individual hydrazines taken from stored drums is shown in Table I.

TABLE I.  
CORRECTION VALUES FOR  $\text{HClO}_4$  TITRATIONS OF HYDRAZINES

Sample No.	$\text{N}_2\text{H}_4$ ml $\text{HClO}_4$ after $\text{Ac}_2\text{O}$	UDMH ml $\text{HClO}_4$ after $\text{Ac}_2\text{O}^1$	MMH ml $\text{HClO}_4$ after $\text{Ac}_2\text{O}$
1	0.20	0.04	0.14
2	0.23	0.03	0.10
3	0.20	0.06	
4	0.21		
5	0.20		
6	0.17		
7	0.05		

<sup>1</sup>ml  $\text{HClO}_4$  less than total alkalinity value. Sample size 0.025 g.

$N_2H_4$  drums 1, 2, and 3 were stored for over five years, while drums 4 and 6 were stored for over two years. Drum 5 was stored only six months. Sample marked 7 was redistilled over a 75 plate distillation column. The UDMH and MMH drums were stored for one year.

The values presented above result from the presence of contaminants in the hydrazines. Aniline and ammonia are the usual titratable contaminants in  $N_2H_4$ . Nitroso-dimethylamine, ammonia, dimethylamine, methylene dimethylhydrazine, and methylhydrazine are the titratable contaminants in UDMH (12). Of these contaminants, aniline and methylhydrazine react similarly to  $N_2H_4$  and are nontitratable after acetylation. Methylene dimethylhydrazine, dimethylamine, ammonia, and nitrosodimethylamine will be acetylated at different rates and could titrate as UDMH. Further work is necessary to establish the rates of reaction for these contaminants. However, their presence in the hydrazines are not significant except for dimethylamine (0.7%) and nitrosodimethylamine (0.3%) in UDMH and ammonia (0.3%) in  $N_2H_4$ .

The contaminants usually present in MMH have not been identified as yet. Until they are, titrations of the individual hydrazines are necessary.

Preliminary observations using other acetylation compounds such as acetyl chloride, benzoyl chloride, maleic anhydride, and benzoic anhydride with UDMH indicate that only the two anhydrides can be used as possible substitutes for acetic anhydride.

#### Effect of Water.

The water content in  $N_2H_4$ -UDMH samples was increased from 0.99% to 7.84% and the samples were analyzed by the procedure presented. Table II shows that increasing the water content does not interfere with the titration.

TABLE II

Percent Water	Sample Weight $\text{N}_2\text{H}_4$ /UDMH before $\text{H}_2\text{O}$ added	Total Weight of Sample	ml $\text{HClO}_4$ Total Alkalinity	ml $\text{HClO}_4$ after Acetylation
0.99	0.8918	0.9007	4.24	1.42
2.31	0.8913	0.9124	4.24	1.42
4.65	0.8903	0.9337	4.24	1.42
7.84	0.8818	0.9570	4.19	1.40

The large excess of acetic anhydride in the  $\text{N}_2\text{H}_4$  reaction reacts with any water to form acetic acid and eliminates the water interference problem encountered in most non-aqueous titration procedures.

#### Potentiometric and Colorimetric Endpoints.

Both the potentiometric and quinaldine red indicator endpoints were evaluated. A Beckman zeromatic pH meter with glass-calomel electrodes was used to determine the stoichiometric endpoints for both the total alkalinity titration and the acetic anhydride titration. Figure III shows a typical curve obtained in these titrations. The quinaldine red endpoint changes from red (actually pink) to colorless approximately 0.01 ml after the greatest change in potential occurs, dictated by the millivolt change plot shown in Figure IV. The indicator and potentiometer endpoints can be improved appreciably by using acetonitrile as the solvent in both titrations. For the plots shown an acetonitrileacetic acid mixture (4:1 V/V) was used.

The endpoints can also be greatly improved by using either dioxane or nitromethane as solvents. However, we have encountered considerable difficulty in obtaining dioxane with a low blank value. The nitromethane is such a good monopropellant that we avoid using it as a solvent except when absolutely necessary.

## PROCEDURE.

### Preparation of Sample.

Pipette 0.4 ml of the  $N_2H_4$ -UDMH mixture into a tared, cooled 50-ml volumetric flask containing 40 ml of acetic acid. Weigh to the nearest 0.1 mg, obtaining the sample weight by difference. Dilute to the 50-ml mark with acetic acid and mix thoroughly.

### Determination of the Total Basicity of the Hydrazine-UDMH Mixture -

Pipette a 5 ml aliquot from the prepared sample into a 50-ml beaker containing 20 ml of acetonitrile. Add 5 drops of quinaldine red indicator and titrate with 0.1 N perchloric acid in acetic acid until the red color of the indicator changes to colorless. Indicate this amount for value "A" in calculations. Titrate a blank of acetonitrile and quinaldine red indicator to the colorless endpoint.

Determination of 1, 1-Dimethylhydrazine - Pipette a 5 ml aliquot of the prepared sample into a 50-ml beaker containing 20 ml of acetonitrile-acetic acid (4:1 V/V). Add 2 ml of acetic anhydride and 5 drops of quinaldine red indicator. Titrate with 0.1 N perchloric acid in acetic until the red color of the indicator changes to colorless. Indicate this amount for value "B" in calculations. Titrate a blank of acetonitrile-acetic acid, acetic anhydride and quinaldine red indicator to the colorless endpoint.

## CALCULATIONS.

$$\% \text{ Hydrazine} = \frac{(A-a) - (B-b) \times N \times M_1}{w} \times 10$$

$$\% \text{ UDMH} = \frac{(B-b) \times N \times M_2}{w} \times 10$$

$$\% H_2O = 100 - (\% \text{ Hydrazine} + \% \text{ UDMH})$$

where:

- A = ml perchloric acid for total basicity.
- B = ml perchloric acid for the anhydride addition.
- a = ml perchloric acid for blank of acetonitrile and indicator.
- b = ml perchloric acid for blank of acetonitrile-acetic acid, acetic anhydride and indicator.
- N = normality of perchloric acid.
- M<sub>1</sub> = milliequivalent weight of hydrazine (0.03205).
- M<sub>2</sub> = milliequivalent weight of UDMH (0.06010).
- w = weight of sample

#### REAGENTS

Acetic acid - glacial reagent grade.

Acetic anhydride - ACS reagent grade.

Acetonitrile - ACS reagent grade.

Quinaldine red - 0.2 percent in acetic acid.

0.1 N perchloric acid in acetic acid - dissolve 8.5 ml of 72% perchloric acid in 1 liter of glacial acetic acid and standardize against primary reagent potassium acid phthalate.

#### RESULTS.

Several samples of N<sub>2</sub>H<sub>4</sub>-UDMH mixtures were prepared and analyzed. Table III presents the theoretical, experimental, and variation from calculated percent obtained for samples ranging from 10% N<sub>2</sub>H<sub>4</sub>-90% UDMH to 90% N<sub>2</sub>H<sub>4</sub>-10% UDMH.

TABLE III

<u>Experimental</u>		<u>Theoretical</u>		<u>Variation in %</u>	
<u>N<sub>2</sub>H<sub>4</sub></u>	<u>UDMH</u>	<u>N<sub>2</sub>H<sub>4</sub></u>	<u>UDMH</u>	<u>N<sub>2</sub>H<sub>4</sub></u>	<u>UDMH</u>
92.18	7.58	92.31	7.33	-0.13	+0.23
96.34	13.03	86.61	13.05	-0.27	-0.02
73.98	25.43	74.31	25.38	-0.33	+0.05
65.59	34.36	65.29	34.67	+0.30	-0.31
55.99	44.35	55.97	44.03	+0.02	+0.32
57.79	42.03	58.18	41.82	-0.39	+0.21
55.12	45.13	55.20	44.80	-0.08	+0.33
55.77	43.80	55.91	44.09	-0.14	-0.29
50.40	48.05	50.30	48.25	+0.10	-0.20
20.65	78.36	20.25	78.80	+0.40	-0.44
12.20	87.63	11.69	88.17	+0.51	-0.54

Standard Deviation N<sub>2</sub>H<sub>4</sub> 0.28

UDMH 0.32

The percent error is greater for samples containing larger UDMH content. This can be attributed to the fact that most of the impurities in the mixtures will react similarly to UDMH. In the above samples of N<sub>2</sub>H<sub>4</sub>/UDMH mixtures, the correction values were not used. The correction values for both the N<sub>2</sub>H<sub>4</sub> and UDMH were 0.04 ml.

Several samples of MMH-UDMH mixtures were prepared and analyzed.

Table IV presents the theoretical, experimental, and variation from calculated value.

TABLE IV

<u>Experimental</u>		<u>Theoretical</u>		<u>Variation in %</u>	
<u>MMH</u>	<u>UDMH</u>	<u>MMH</u>	<u>UDMH</u>	<u>MMH</u>	<u>UDMH</u>
41.56	56.54	41.73	57.04	-0.17	-0.50
50.42	49.24	51.06	48.19	-0.64	+1.05
40.07	55.88	40.47	56.06	-0.40	-0.18
47.01	49.84	47.64	51.19	-0.63	-0.35
62.97	34.47	63.12	34.92	-0.15	-0.45
63.57	36.07	63.89	35.34	-0.32	+0.73

Standard Deviation MMH 0.43  
UDMH 0.61

When the correction values discussed in Table I are made to MMH/UDMH samples, the following results are obtained:

TABLE V

<u>Experimental</u>		<u>Theoretical</u>		<u>Variations in %</u>	
<u>MMH</u>	<u>UDMH</u>	<u>MMH</u>	<u>UDMH</u>	<u>MMH</u>	<u>UDMH</u>
41.85	57.10	41.73	57.04	+0.12	+0.06
50.78	48.06	51.06	48.19	-0.18	0.13

We encounter difficulty in the analysis of the hydrazine mixtures when the correction values for the individual hydrazines are unknown. The correction values are applied by equating the value for the individual hydrazine proportional to the total concentration hydrazines in the mixture. The UDMH correction value (resulting from MMH type contaminants) is added to the milliliters acid required for either the  $N_2H_4$  or MMH titration. The  $N_2H_4$  or MMH correction value (resulting from UDMH type contaminants) is subtracted from the milliliters acid required for the UDMH titration.

When the correction values become excessive, the results will be unsatisfactory. Values this high are not usually encountered in stored hydrazines. We presume that the high-value MMH samples reported previously were contaminated accidentally and are not representative correction values normally encountered.

#### SUMMARY.

Detailed procedures for potentiometric or colorimetric titration in acetonitrile/acetic acid solvent with perchloric acid titrant are presented. Mixtures of either UDMH/ $N_2H_4$  or UDMH/MMH can be analyzed in this system. The procedural basis is the difference in reaction rates for acetylating hydrazines with acetic anhydride.



Analysis of  $\text{N}_2\text{H}_4$ /UDMH mixtures gave standard deviations of 0.28 for  $\text{N}_2\text{H}_4$  and 0.32 for UDMH. Standard deviations found for analysis of MMH/UDMH were 0.43 for MMH and 0.61 for UDMH. The standard deviation decreases when correction values are applied.

#### RECOMMENDATION.

The acetylation procedure for analysis of  $\text{N}_2\text{H}_4$ /UDMH is recommended to supersede the aldehyde procedure currently required in Specification MIL-P-27402.




Figure 1. Effect of Time and Temperature on the Reaction Rate of UDMH with Acetic Anhydride.




Figure II. The Effect of Acetic Anhydride on Hydrazines




Figure III. Typical Titration Curve for Both Titrations.

## BIBLIOGRAPHY

1. Cain, E. F. C. and Stevens, M. R., paper from 2nd International Biannual Symposium on Gas Chromatography (June 1959).
2. Clark, C. C., "Hydrazine," p.43, Copyright (1953) by Mathieson Chemical Corporation.
3. Fritz, J. S. and Fulda, M. O., Anal. Chemistry 25, 1837 (1953).
4. Fritz, J. S. and Schenk, G. H., Anal. Chem. 31, 1808 (1959).
5. Mackey, M. R., unpublished paper, 28 May 1962.
6. Malone, H. E., Anal. Chemistry, 33, 574 (1961).
7. Sadtler, Samuel P. and Son, Inc., Standard Spectra, No. 19955, Aldrich Chemical Co.; 11889, Iowa State College; 11890, Iowa State College.
8. Streuli, C. A., Anal. Chemistry 30, 997 (1958).
9. Weed, A. F., Analytical Method CA-106 by Food Machinery and Chemical Corporation, May (1960).
10. Wimer, D. C., Anal. Chemistry 30, 77 (1958).
11. Wimer, D. C., *ibid*, 30, 2060 (1958).
12. Zarembo, J. D., Analytical Method CA-196 by Food Machinery and Chemical Corporation, February (1962).

## DISTRIBUTION

Rocket Research Laboratories  
Attn: DGPCS, Mr. Hugh Malone  
6593 Test Group (Dev)  
Edwards, Calif (50)

Commander  
AFSC  
Foreign Technology Div  
Attn: FTD (TD-E3b)  
Wright-Patterson AFB, Ohio (1)

Director, Air University Library  
Maxwell AFB, Alabama (2)

Air Force Systems Command  
Attn: SCR-2  
Andrews AFB  
Wash 25, DC (2)

Bureau of Ordnance  
Department of the Navy  
Wash 25, DC (1)

U. S. Naval Research Laboratory  
Attn: Dr. H. W. Carhart  
Wash 25, DC (1)

LPIA  
Applied Physics Laboratory  
Johns Hopkins University  
8621 Georgia Ave  
Silver Spring, Maryland (10)

Department of the Air Force  
Hq USAF, DCS/D  
Attn: AFDRT/AS  
Wash 25, DC (1)

Bureau of Aeronautics  
Department of the Navy  
Attn: PP-731  
Wash 25, DC (2)

MAAMA (MAOQ/Mr. Alvin Traynor)  
Olmsted AFB, Pa. (10)

Office of Naval Research  
Attn: Mr. B. Hornstein  
Washington 25, DC (3)

Office of Chief of Ordnance  
Department of the Army  
Attn: Dr. C. M. Hudson (ORDTU)  
Wash 25, DC (3)

Commanding General  
Army Chemical Center  
Attn: Dr. E. H. Krackow  
Edgewood, Maryland (2)

Picatinny Arsenal  
Attn: Mr. L. H. Rosen  
Dr. J. D. Clark  
Dover, New Jersey (2)

AFMTC (Mr. Normile)  
Patrick AFB, Fla (1)

U. S. Naval Ordnance Test Station  
Attn: Mr. E. M. Bens  
China Lake, Calif (2)

Commander  
Research and Test Squadron  
Attn: Mr. Quirmback  
Mr. M. D. Weber  
Holloman AFB, New Mexico (2)

Naval Powder Factory  
Attn: Mr. W. O. Brimijoin  
Indianhead, Maryland (2)

New York University  
Department of Chemical Engineering  
Attn: Dr. P. F. Winternitz  
New York 53, NY (1)

Distribution p2 of 3

Westvaco Chlor-Alkali Division  
Food Machinery  
Attn: Dr. Paul Derr  
Dr. Louis Diamond  
Princeton, NJ (2)

Commanding General  
Engineering Branch  
Technical and Engineering Div, O. M. L.  
Attn: Mr. D. I. Graham, Jr.  
Redstone Arsenal  
Huntsville, Alabama (1)

Commanding General  
Mechanical Branch  
Technical and Engineering Div, O. M. L.  
Attn: Mr. P. E. Redding  
Redstone Arsenal  
Huntsville, Alabama (2)

Aerojet-General Corporation  
Attn: Dr. D. L. Armstrong  
Mr. E. M. Wilson  
Azusa, Calif (2)

Bell Aircraft Corporation  
Attn: Dr. S. A. Long  
Mr. Henry Heubusch  
PO Box 1  
Buffalo, NY (3)

Phillips Petroleum Company  
Attn: Mr. Homer Fix  
Bartlesville, Okla (2)

Commercial Solvents Co.  
Research and Development Div  
Attn: Dr. R. S. Egly  
Dr. J. A. Riddick  
Terre Haute, Indiana (2)

Olin-Mathieson Chemical Corporation  
Attn: Mr. D. Griffin  
New Haven, Conn (2)

Olin-Mathieson Chemical Corporation  
Attn: Dr. B. Hill  
Baltimore, Maryland (2)

General Chemical Division  
Allied Chemical and Dye Corporation  
Attn: Mr. T. A. Wallace  
40 Rector St  
New York, NY (1)

Southern Research Institute  
Attn: Dr. W. J. Barrett  
917 S Twentieth Ave  
Birmingham, Ala (1)

California Institute of Technology  
Jet Propulsion Laboratory  
2800 Oak Grove Dr  
Pasadena, Calif (1)

Commander  
Air Proving Grounds, Chemical Section  
Climatic Project Air Corps  
Attn: Mr. E. G. Taylor  
Eglin AFB, Fla (1)

General Electric Company  
Rocket Propulsion Units  
Attn: Mr. E. St Clair Gantz  
Building 300  
Cincinnati 15, Ohio (1)

Aerojet-General Corporation  
Attn: Dr. W. R. Fish  
Sacramento, Calif (1)

Aerojet-General Corporation  
Attn: Mr. Harry York  
Sacramento, Calif (1)

Rohm and Haas Co.  
Redstone Div  
Attn: Dr. Keith A. Booman  
Huntsville, Ala (1)

Rocketdyne  
Attn: Mr. E. F. C. Cain, Code  
D/591-356  
6633 Canoga Ave  
Canoga Park, Calif (1)

Distribution p3 of 3

New York University  
College of Engineering, Bldg 3  
Attn: Dr. Arrigo A. Carotti  
New York 53, NY (1)

Naval Research Lab, Code 6180  
Attn: Dr. Robert Hazlett  
Wash 25, DC (1)

Thiokol Chemical Corporation  
Reaction Motors Div  
Attn: Mr. Fred G. Hoffman  
Denville, NY (1)

George C. Marshall  
Space Flight Center  
NASA M-S and M-MA  
Attn: Mr. John Nunnelley  
Huntsville, Ala (1)

Lewis Research Center  
Propellant Chemistry Div  
Analytical Chemistry Sec  
Attn: Mr. Adolph Spakowski  
Cleveland, Ohio (1)

General Chemical Co.  
Research Laboratory  
Attn: Mr. A. D. Turissini  
Morristown, NJ (1)

Army Rocket and Guided Missile Agency  
Research Laboratory, OMLD  
Attn: Dr. Walter Wharton,  
Code ORDXR-ORC  
Huntsville, Ala (1)

Allied Chemical Corporation  
Nitrogen Div  
Attn: Dr. R. E. Ogden  
Hopewell, Virginia (1)

Air Products and Chemicals Inc  
Attn: Mr. Willard Ent  
Allentown, Pa. (1)

Texaco Research Center  
Attn: Dr. H. E. Vermillion  
Beacon, NY (1)  
AFSC (STLO)  
Lewis Research Center (NASA)  
Cleveland 35, Ohio (1)

United Technology Corporation  
Attn: Mr. Paul Marti  
PO Box 358  
Sunnyvale, Calif (1)

The Martin Company (C-124)  
Research Dept  
Attn: Dr. James A. Bowman  
PO Box 1176  
Denver 1, Colo (1)

Stanford Research Institute  
Poulter Laboratories  
Attn: Dr. Eugene Burns  
Menlo Park, Calif (1)

Office of Naval Research  
Department of the Navy  
Attn: Mr. Thomas Dowd  
Wash 25, DC (1)

AFLC (MCRQD, Mr. Russell Dye)  
Wright-Patterson AFB, Ohio (1)

AFSC (MCQDT/Mr. Stephen Dzuroff)  
Wright-Patterson AFB, Ohio (1)

Armed Service Tech Information  
Agency  
Arlington Hall Station  
Attn: TIPCR  
Arlington 12, Virginia (10)

Air Force Flight Test Center  
Attn: FTOOT  
Edwards AFB, Calif (10)

Rocket Research Laboratories  
6593d Test Group (Dev), (DGLP)  
Edwards, Calif (10)

DCAS (DCLMT/Mrs. Densmore)  
AF Unit Post Office  
Los Angeles 45, Calif (1)

BSD (BSREQ, Dr. Paul Leatherman)  
Norton AFB, Calif (1)

Robert A. Biggers (DGPCS)  
Rocket Research Laboratories  
6593 Test Group (Dev)  
Edwards, Calif (1)



<p>Rocket Research Laboratories, Edwards, Calif. Rpt No. SSD-TDR-62-203. AN ACETYLATION METHOD FOR DETERMINING MIXTURES OF EITHER HYDRAZINE-1, 1-DIMETHYLHYDRA- ZINE OR MONOMETHYL-HYDRAZINE-1, 1-DIMETHYLHYDRAZINE. 14 p. incl illus. Dec 62.</p> <p>Unclassified Report</p> <p>A method is presented for determining mixtures of either hydrazine-1, 1-dimethylhydrazine (N<sub>2</sub>H<sub>4</sub>/UDMH) or monomethylhydrazine-1, 1-dimethylhydrazine (MMH/UDMH). This method is based on the difference in reaction rates for acetylating hydrazines. Both N<sub>2</sub>H<sub>4</sub> and MMH react immediately while the acetylation of UDMH proceeds slowly. The reactions are conducted in</p>	<p>1. Acetylation</p> <p>I. AFSC Project No. 3148, Task No. 30197</p> <p>II. Malone, Hugh E., Biggers, R. A.</p> <p>III. In ASTIA collection</p>
---	--

<p>Rocket Research Laboratories, Edwards, Calif. Rpt No. SSD-TDR-62-203. AN ACETYLATION METHOD FOR DETERMINING MIXTURES OF EITHER HYDRAZINE-1, 1-DIMETHYLHYDRA- ZINE OR MONOMETHYL-HYDRAZINE-1, 1-DIMETHYLHYDRAZINE. 14 p. incl illus. Dec 62.</p> <p>Unclassified Report</p> <p>A method is presented for determining mixtures of either hydrazine-1, 1-dimethylhydrazine (N<sub>2</sub>H<sub>4</sub>/UDMH) or monomethylhydrazine-1, 1-dimethylhydrazine (MMH/UDMH). This method is based on the difference in reaction rates for acetylating hydrazines. Both N<sub>2</sub>H<sub>4</sub> and MMH react immediately while the acetylation of UDMH proceeds slowly. The reactions are conducted in</p>	<p>1. Acetylation</p> <p>I. AFSC Project No. 3148, Task No. 30197</p> <p>II. Malone, Hugh E., Biggers, R. A.</p> <p>III. In ASTIA collection</p>
---	--

<p>acetonitrile-acetic acid medium. Two titrations are required; one measures the total basicity of the hydrazine mixtures; the other titration meas- ures UDMH after the N<sub>2</sub>H<sub>4</sub> or MMH has reacted with Ac<sub>2</sub>O.</p>	
---	--

<p>acetonitrile-acetic acid medium. Two titrations are required; one measures the total basicity of the hydrazine mixtures; the other titration meas- ures UDMH after the N<sub>2</sub>H<sub>4</sub> or MMH has reacted with Ac<sub>2</sub>O.</p>	
---	--